



Mount Sinai - Elmhurst
Emergency Critical Care

RACC PRE-SHIFT CHECKLIST

AIRWAY

1. **Box check:**
 - a. MAC 3 and 4 check: Laryngoscope bulb functional, screwed in tightly
 - b. ETT 8.0 balloon intact, stylet in place, properly shaped. 7.5 ETT open, balloon intact.
 - c. Colorimetric EtCO₂ package sealed
 - d. Oral airway open
 - e. This is now your box until you use it. Reseal it with tape, label with your name and date
2. **Video laryngoscope cart check:**
 - a. CMAC size 3 and 4 blades
 - b. Hyperangulated D blade, rigid stylet available
 - c. Failed airway kit check:
 - i. LMA 3.5 and 4.5, packaging open, syringe attached
 - ii. Scalpel No.11
 - iii. Portex 6.0
 - d. Bougie x4

HEADWALL

1. Nasal cannula
2. NRB
3. Face shields
4. PEEP valve
5. Oral airway
6. BVM taken out of bag, connected to O₂, hung back on regulator
7. Suction available
8. Sidestream EtCO₂ connected, functional
9. Zoll defibrillator:
 - a. pads connected
 - b. EtCO₂ adapter
10. Ventilator:
 - a. Tubing/filter connected

RACC INVENTORY CHECK

When stock is less than 2/3rds inventory, please contact appropriate department and document on inventory sheet. If department is not available or item not obtained, signout to oncoming team.

1. Airway cart inventory
2. Critical resuscitation inventory
3. Respiratory equipment inventory
 - a. BiPAP masks
 - b. Nasal cannula
 - c. NRB
 - d. Nebulizers, handheld and mask
 - e. Venturi masks
 - f. PEEP valves
 - g. Bag valve masks