



Mount Sinai - Elmhurst  
Emergency Critical Care

# TRAUMA PRE-SHIFT CHECKLIST

## AIRWAY

1. **Box check:**
  - a. MAC 3 and 4 check: Laryngoscope bulb functional, screwed in tightly
  - b. ETT 8.0 balloon intact, stylet in place, properly shaped. 7.5 ETT open, balloon intact.
  - c. Colorimetric EtCO<sub>2</sub> package sealed
  - d. Oral airway open
  - e. This is now your box until you use it. Reseal it with tape, label with your name and date
2. **Video laryngoscope cart check:**
  - a. CMAC size 3 and 4 blades
  - b. Hyperangulated D blade, rigid stylet available
  - c. Failed airway kit check:
    - i. LMA 3.5 and 4.5, packaging open, syringe attached
    - ii. Scalpel No.11
    - iii. Portex 6.0
  - d. Bougie x4

## HEADWALL

1. Nasal cannula
2. NRB
3. Face shields
4. PEEP valve
5. Oral airway
6. BVM taken out of bag, connected to O<sub>2</sub>, hung back on regulator
7. Suction available
8. Sidestream EtCO<sub>2</sub> connected, functional
9. Zoll defibrillator:
  - a. pads connected
  - b. EtCO<sub>2</sub> adapter
10. Ventilator:
  - a. Tubing/filter connected

## TRAUMA CART

1. Chest tubes
2. Thoracostomy tray
3. Gauze 4x4s, betadine
4. Code meds
5. Kelley clamps for pelvic binder

## TRAUMA DIVIDER

1. LMA 3.5 and 4.5, 12cc syringe
2. Suture kit x2
3. Large skin stapler x2

## RESUSCITATION EQUIPMENT

1. Critical resuscitation equipment inventory
2. Trauma resuscitation equipment inventory
3. Respiratory equipment inventory
  - a. BiPAP masks
  - b. Nasal canulas
  - c. NRB
  - d. Nebulizers, handheld and mask
  - e. Venturi masks
  - f. PEEP valves
  - g. Bag valve masks