

START OF RESUS CHECKLIST

Airway Box: check only 1 box

- ETT 8.0 and ETT 7.5 balloons intact, stylet in place, properly shaped
- Oral airway adjuncts present
- ETT securing device present
- Bougie present
- Reseal the box with tape, label with your name and date
- To replace airway boxes page #4728 or call 4-PEP. 3 Airway boxes should be present.

General Airway Prep

- Check BIPAP masks: small, medium, large
- Check that PEEP valves are stocked
- Check that non-rebreather masks, nasal cannula, and venti masks with albuterol attachment are stocked
- Check that LMA are present
- Check that perc cric kits are stocked
- Check glidescope is present and works
- Check glidescope stylets are stocked
- Check glidescope disposable covers
- Check #11 scalpels & keep one on hand

Resus Bays

- Check that BVM is present and connected to oxygen source
- Check wall oxygen and ensure Y-adaptor is attached.
- Set-up and check suction with Yankauer catheter and tubing x2

- Set-up and check ETCO2 monitor, tubing, and water trap
- Ensure vent is present and tubing with filter is connected

Access

- Check vessel finder
- Check that long IV catheters in 20, 18, and 16 gauge sizes are stocked
- Check that triple lumen central venous access kits are stocked
- Check that percutaneous introducer kits are stocked
- Check transvenous pacer set-up is stocked
- Check chlorhexidine is stocked
- Check IO gun, needles.

General Resus Prep

- General imaging US are present, working and clean
- Check that defibrillators are present, working and clean
- Check that chest tubes are stocked
- Check that pigtail catheter kits are stocked
- Check that pleuravac set-up is stocked
- Check that chest trauma tray is stocked

If any equipment or supplies are missing, notify charge RN at 24015